

Task 3 Report: Feature Extraction and Comparison from 3D Knee CT Scan

Objective

To design and implement a complete pipeline for 3D feature extraction and comparison from knee CT scans, segmented into anatomical regions: **Tibia**, **Femur**, and **Background**. The features are extracted using a **3D-inflated version of a pre-trained 2D DenseNet121**, and their similarities are evaluated using **cosine similarity**.

Pipeline Components

1. Segmentation-Based Splitting

- A **3D segmentation mask** is loaded and used to isolate three anatomical regions:
 - **Tibia (green)** – label 1
 - **Femur (red)** – label 2
 - **Background** – all other voxels
- The CT volume and the mask are aligned to extract per-region volumes for analysis.

2. 2D to 3D Model Conversion

- A pretrained **2D DenseNet121** is loaded from `torchvision.models`.
- All **Conv2D** layers are **inflated into Conv3D**:
 - The 2D weights (`out_ch, in_ch, H, W`) are repeated along a new depth axis \rightarrow (`out_ch, in_ch, D, H, W`)
 - The weights are normalized by dividing by `D` to maintain activation scale.

3. Feature Extraction

- Each **3D region volume** (Tibia, Femur, Background) is passed through the **3D DenseNet121**.
- Feature maps are extracted from:
 - Last convolutional layer
 - Third-last convolutional layer
 - Fifth-last convolutional layer
- **Global Average Pooling (GAP)** is applied to convert each feature map to a fixed-length feature vector.

4. Feature Comparison

- Pairwise **cosine similarity** is computed between feature vectors:
 - Tibia \leftrightarrow Femur
 - Tibia \leftrightarrow Background
 - Femur \leftrightarrow Background
- Similarity is computed **separately for each of the 3 feature levels** (last, third-last, fifth-last).

5. Result Organization

- Cosine similarity scores are **saved into a CSV file**.
 - **Rows:** Each image (or region pair)
 - **Columns:** Similarity for each of the 3 layers
- Example format:

Pair	Fifth-Las st	Third-Las st	Las t
Tibia-Femur	0.78	0.84	0.88
Tibia-Background	0.32	0.36	0.41
Femur-Background	0.28	0.34	0.39

Improvements with More Time

1. Fine-Tuning the 3D CNN on Medical Data

- **Current Limitation:** The DenseNet121 model is pre-trained on natural images (ImageNet), which are quite different from medical scans.
- **Improvement:** Fine-tune the inflated 3D DenseNet on a medical imaging dataset (e.g., **MedicalNet**, **MRNet**, or custom annotated knee CTs) to

improve feature relevance and accuracy.

2. Use of True 3D Pretrained Models

- Instead of inflating 2D models, leverage pretrained 3D models such as:
 - **MedicalNet** (pretrained on multiple medical datasets)
 - **MONAI**-based networks
 - **3D ResNet**, **3D DenseNet** from [torchvision](#) or [video_classification](#)

3. Segmentation Model Integration

- Currently, masks are provided.
- Add a step to **auto-segment the CT scans** using:
 - 3D U-Net
 - nnU-Net

4. Augmentation and Preprocessing Enhancements

- Add 3D data augmentation techniques to increase generalization:
 - Random crop

- Rotation
- Elastic deformation
- Intensity scaling
- Normalize intensities using medical-specific heuristics (e.g., clipping Hounsfield Units).